

THE KEYWING

An award-winning key turner designed for those with arthritis or weak hands



This award winning design offers a simple solution to the problem of using Yale Keys which have a small surface area making them difficult to grip and turn. The Keywing clips over the key to provide a much wider finger contact area which improves grip and leverage when turning. Available in a range of 3 colours allowing the user to colour-code their keys and find the right one first time, every time.

HOW TO USE:

- 1. Pull clips to release tray
- 2. Peel off the film and stick your key to the pad
- 3. Reinsert the key and tray until it clicks into place

FEATURE	BENEFIT		
Increases size of key	Easier to grip and provides leverage for turning		
Bright colours available	Easy to find		
Tactile bumps	Stop fingers slipping		
Made from Nylon	Warm to touch even in extreme cold		

RECOMMENDED FOR:

KAN			
Weak Grip	Poor Hand Control	Visual Impairment	Tremor

SIZES & CODES:

Single: 170 x 95 x 10 mm 13g (inc packaging)

PAP-KW2BL – Blue PAP-KW2GN – Green PAP-KW2PK – Pink **Pack of 3**: 170 x 95 x 10 mm 28g (inc packaging) PAP-KW2X3 – Blue/Green/Pink

Nov-2022