

# PRIMA OUTDOOR GRAB BAR

A WELL SITUATED GRAB BAR CAN MAKE ALL THE DIFFERENCE TO THE SAFETY AND CONFIDENCE OF USERS OUTSIDE THEIR HOME.

## KEY FEATURES

- Soft ribs help to improve grip, even when wet
- Strong plastic bar which is warm to the touch, unlike metal
- Concealed aluminium core for a high maximum user weight

The Prima Outdoor Grab Bar is the ideal rail for outside the front or back door, helping users get in and out of their home.

The subtle colour is a lot more discreet than white metal rails and can blend into its surrounding when positioned on brick work.

It is fully weather resistant, and features our unique Prima rubber grips which help to prevent the user's hand slipping, even when wet. It is always warm to the touch, and has an oval cross-section which ergonomically matches the human grip.

Grab bars can be subject to some hefty pulls and sudden tugs, so they are made from robust materials and have strong hidden fixing points. The Prima Outdoor Grab Bar has a high maximum load limit of 160kg / 25 stone, due to a concealed aluminium core running through the middle of the bar.

It is essential that grab bars are fitted competently, according to the instructions, or seek professional advice from a qualified tradesperson.



6537

Single Retail Pack

6538

Bulk Pack (order in 20's)



**160KG / 25 STONE / 350LBS**

