

△ HelpSoq ▽



Brand

Sub-brand

GS1 article code (GTIN)

Product classification (GS1/GTIN)

Article description (GS1/GTIN)

Function

Box Contents

Controls

Accessories

Box Dimensions (Exterior)

Device Dimensions

Total weight

Device weight

Materials

Electrical Requirements (Input)

Engine

Designed by

Produced and assembled in

CE Directives

TÜV

Patents

Manufacturer

Contact info Email

Contact info Phone

Not suitable for

HelpSoq On

HS01

8720865049038

10008118 - Supporting Component of a Medical Device

Donning aid for TEK sizes M/L/XL

Electric aid for donning Therapeutic Elastic Stockings (TEK)

HelpSoq On device, Power supply, Power cord, Manual, Instruction sheet, Booklet: "I want compression stockings"

Slide switch for Off, Pause, and On. Red Reverse wheel rotation Button

Ring (two-piece, diameter 14 cm; circumference 45 cm), enclosed

Length x width x height : 26 x 26 x 31 cm

height: 23 cm; diameter: 22 cm; circumference : 70 cm

5.2 kg. Includes device, power supply, cord, documentation, and packaging

2,5 kg

Plastics (no Latex), stainless steel, magnets, and printed circuit board (PCB) with electronic components

AC 100-240V, 2.5A, 50-60Hz mains power supply (enclosed)

DC 24.0V - 4.5A

HelpSoq B.V.

The Netherlands

Complies with the relevant provisions of Directives: LVD 2014/35/EU; EMC

2014/30/EU; RoHS 2011/65 EU

In progress

Dutch patent 2025439, Int'l patent pending PCT/NL2021/050174

HelpSoq B.V., Hondiusstraat 28 C, 6827 DE, Arnhem (NL)

info@helpsoq.com

+31 85 400 7744

TEK stockings with a cB value of 21 cm or smaller (commonly known as size xS, S, Small, or 'l') can be damaged by extreme stretching;

also not suitable for normal stockings or TEK stockings Class 0 and 1