

Let's Move Again

Outdoor rollator

Green - Black Art. No. LMA601 EAN: 07350036885298
Grey - Black Art. No. LMA701 EAN: 07350036885304



Let's Move Again, stylish and innovative outdoor rollator from Trust Care. The rollator has a sturdy lightweight aluminum frame and whole rollator weighs only 5.4 kg. With such low weight it is easy and pleasurable to use the rollator in any conditions. Of course it has an effective brake and parking brake that are easy to operate.

A fabric bag is also available for Let's Move Again. The rollator is easy to fold and can stand by itself when folded. Folded width is only 20 cm, such feature makes it easy to store or take it with you wherever you go. Let's Move Again has ergonomically designed rubber handles, adjustable from 78 to 93 cm in order to fit different users.

Let's Move Again is tested and approved according to ISO specifications for 130 kg user weight.



Technical specifications:

Tested according to ISO 11199-2:2005

Max. user weight	130 kg
Width	58 cm
Length	65 cm
Weight	5.4 kg
Handle height is adjustable from	78 – 93 cm
Seat height	60 cm
Wheel diameter	20 cm
Max. bag load	12 kg

